



Hycult Scope



Complement markers of renal disease and injury

Complement is a major mediator system in pathogenesis of various kidney diseases. The presence and localization of complement components in glomerulus and/or the tubulo-interstitial area provides diagnostic tools for several human renal diseases (Table).

Determination of complement deposition also serves a valuable supplementary role in diagnosis of antibody mediated rejection in kidney transplantation. Particularly important has been the inclusion of C4d staining in the Banff scoring criteria for assessment of renal transplant pathology, in 2003.

In addition, the application of anti-complement antibodies is of utmost importance for studying experimental models of kidney disease in rat and mouse. A specific example is the mouse model where antibodies directed against components of the globular basement membrane (GBM) cause renal injury (Figure 1b). Deposition of C1q along the GBM demonstrates an activated innate immune response.

Furthermore, localization studies of complement proteins provide insight into the underlying pathophysiology of acute kidney injury. Staining for complement proteins has therefore become a routine part of renal biopsy analysis. However, early diagnosis of renal tissue damage requires non-invasive biomarkers, like NGAL (neutrophil-gelatinase-associated lipocalin), nitrotyrosine and FABPs (fatty acid binding proteins) that can easily and timely be measured in serum or urinary samples.

DETERMINANTS OF HUMAN RENAL DISEASE AND INJURY

| Disease | Complement deposition / Biomarkers |
|---|--|
| Poststreptococcal glomerulonephritis | C3 ^a |
| Membranous glomerulonephritis | C3 ^a |
| Membranoproliferative glomerulonephritis I | C1q ^a , C4 ^a , C3 ^a |
| Membranoproliferative glomerulonephritis II | C3 ^a |
| IgA nephropathy | C3 ^a , {C4d ^a , MBL ^a , MASP ^a , L-ficolin ^a } [†] |
| Lupus nephritis | C1q ^a , C4 ^a , C3 ^a |
| Antibody-mediated rejection | C4d ^b |
| Acute kidney injury (AKI) | C3d ^c , C8 ^d , MBL ^{b,c} / NGAL, Nitrotyrosine, FABP |

a = in glomeruli; b = in peritubular capillaries; c = in renal tubules; d = in tubulointerstitium; † = in subpopulation of patients





Figure 1 A) Glomerulus. B) Deposition of mouse C1q in an anti-GBM model of glomerulonephritis. Glomerulonephritis was induced in mice by injection of Rabbit anti-Mouse Glomerular Basement Membrane (GBM) IgG. Activation of complement is shown by C1q deposition along the GBM (Cat.# HM1096).

References:


Coca et al., *Kidney Intern.* (2007)
 Molitoris et al., *Nat Clin Pract Nephrol* (2008), 4:154-165
 Portilla et al., *Kidney Intern* (2008), 73: 465-472
 Roos et al., *Nephrol Dial Transplant* (2007), 22: 3370-3377
 Thurman et al., *Kidney Intern.* (2005), 67: 524-530

Thurman et al., *Kidney Intern.* (2008), 73: 379-381
 Trouw et al., *Mol Immunol* (2003), 40: 125-134
 de Vries et al., *Am J Pathol* (2004), 165: 1677-1688
 de Vries et al., *Transplantation* (2003), 75: 375-382
 Yamamoto et al., *J Am Soc Nephrol* (2007), 18: 2894-2902
 Zager et al., *Am J Renal Physiol* (2006), 291: 546-556

Kidney diseases

In this Issue

| | |
|--|---|
| Complement markers of renal disease and injury | 1 |
| Urinary FABPs as biomarkers in renal disease | 2 |
| NGAL: biomarker for renal injury and neoplasia | 2 |
| The role of complement in ANCA related necrotizing small vessel vasculitis | 3 |
| Soluble MAdCAM-1, an inflammatory biomarker | 4 |



Urinary FABPs as biomarkers in renal disease

J.M.Hofstra MD, Dept. of Nephrology, Radboud University Nijmegen Medical Center, Nijmegen, The Netherlands

Renal disease is a common health problem. Progression of renal injury leads to end stage renal failure, which requires renal replacement therapy (hemodialysis, kidney transplantation). The recognition that tubular injury is important in progressive renal injury has stimulated the interest in urinary markers of tubular cell injury. Several tubular markers, such as β -2-microglobulin, α -1-microglobulin, Alkaline Phosphatase or β -NAG can predict prognosis in patients with renal diseases. The above-mentioned proteins are rather non-specific markers of proximal tubular cell injury. Little is known on the role of distal tubular injury in progressive renal disease. In addition, more specific or pathogenetically relevant markers may offer advantages in evaluating various patient groups.

Fatty-acid binding proteins (FABPs) are intracellular carrier proteins. Two types of FABPs are localized in human renal tubular cells. L-FABP is found in the cytoplasm of proximal tubules, whereas H-FABP is localized in the distal tubules. An increased urinary excretion of FABPs may simply result from release by structurally damaged tubular cells. However, it is suggested that tubular L-FABP expression is upregulated by hypoxia and increased excretion may thus occur before the occurrence of the actual structural damage, making it an early biomarker. In our patients with membranous nephropathy both L-FABP and H-FABP predicted prognosis with rather high sensitivity and specificity. Further investigation in differences between proximal and distal tubular markers like L-FABP and H-FABP may bring new insights to the pathogenesis and process of tubulo-interstitial damage.

HOT NEWS IN KIDNEY RESEARCH:

Akirin

Akirin (Japanese for "making things clear") plays an important role in triggering the innate immune response of the fruit fly *Drosophila*, mice and even humans. Suppressing Akirin production in immune cells of the fly leads to a significantly enhanced susceptibility to bacterial infections.

Goto et al Nature Immunology 2008, 9: 97-104

KIM-1

Kidney injury molecule-1 (KIM-1) staining identifies sensitively and specifically proximal tubular injury and correlates with degree of renal dysfunction.

Zhang et al. Kidney Int 2008, 73: 608-14

Get "stoned" from fructose


High fructose intake is associated with increased risk of kidney stones. Although patients with nephrolithiasis require increased fluid consumption, the authors advise them to avoid intake of fructose-rich drinks like sugar sweetened soft drinks, diet soft drinks, and fruit juice.

Taylor et al. Kidney Int. 2008, 73: 207-12

NGAL: biomarker for renal injury and neoplasia

A sensitive, non-invasive tool for monitoring epithelial injury and carcinogenesis (Cat. # HK330)

- human neutrophil gelatinase-associated lipocalin (NGAL)
- induced in most tissues exposed to microorganisms, and in epithelial cells during inflammation
- antimicrobial properties and plays a role in regulating inflammation and cellular growth
- detectable in several body fluids upon epithelial injury
- important and informative biomarker of acute ischemic injury
- biomarker for early detection of renal failure
- NGAL activates nephron formation in the embryonic kidney, is rapidly and massively induced in renal failure and possesses kidney-protective activity
- blood, urine, and kidney NGAL levels are the real-time indicators of active kidney damage, rather than one of many markers of functional nephron number
- biomarker in ovarian cancer, colorectal cancer, adenocarcinomas, urothelial carcinomas, ischemic cerebrovascular disease



Kidney damage

ASSAYS

| SPECIFICITY | CAT. # |
|--|---------------|
| H-FABP, Human | HK401 |
| H-FABP, Mouse, Rat | Unique HK403 |
| NGAL, Human | HK330 |
| Nitrotyrosine | Unique HK501 |
| Lectin NPN early apoptosis detection kit | Unique HIT303 |
| α -Defensin 1-3, Human | Unique HK317 |
| Annexin V-FITC | HIT304 |
| Calprotectin, Human | HK325 |
| Elafin/SKALP, Human | Unique HK318 |
| Elastase, Human | HK319 |
| LL-37, Human | Unique HK321 |

MONOCLONAL ANTIBODIES

| SPECIFICITY | CAT. # |
|--------------------------|--------|
| H-FABP, Human, mAb 66E2 | HM2016 |
| H-FABP, Human, mAb 67D3 | HM2018 |
| L-FABP, Human, mAb K5A6 | HM2051 |
| L-FABP, Human, mAb L2B10 | HM2049 |
| NGAL, Human, mAb 697 | HM2193 |
| Nitrotyrosine, mAb HM11 | HM5001 |




Antibodies for Kidney Research

MONOCLONAL ANTIBODIES

| SPECIFICITY | CAT. # |
|---------------------------------------|--------|
| Activated C3, Human, mAb bH6 | HM2168 |
| Activated Protein C, Human, mAb PC107 | HM2151 |
| α-Defensin 1-3, Human, mAb D21 | HM2058 |
| α-V-beta3 Integrin, Human, mAb BV3 | HM2034 |
| Arginase type-1, Human, mAb 6G3 | HM2162 |
| Arginase type-1, Human, mAb 9C5 | HM2163 |
| Barmotin/7H6 antigen, Human, mAb 7H6 | HM2102 |
| β Catenin, Human, mAb 9F2 | HM2112 |
| β Integrin, Human, mAb BV7 | HM2033 |
| C1q, Mouse, mAb 7H8 | HM1044 |
| C1q, Mouse, mAb JL-1 | HM1096 |
| C3, Mouse, mAb 11H9 | HM1045 |
| C3/C3a, C-terminus, Human, mAb 2898 | HM2075 |
| C3/C3a, Human, mAb 474 | HM2073 |
| C3/C3b, Human, mAb 755 | HM2072 |
| C3/C3b, Rat, mAb 2B10B9B2 | HM3031 |
| C3a, Mouse mAb 3/11 | HM1072 |
| C3a/C3a des-Arg, Human, mAb 2991 | HM2074 |
| C3aR, Human, mAb 17 | HM2195 |
| C3aR, Rat, mAb 74 | HM3028 |
| C3b/iC3b/C3c, Mouse mAb 2/11 | HM1065 |
| C3b/iC3b/C3c, Mouse, mAb 3/26 | HM1078 |
| C3c, Human, mAb 4 | HM2200 |
| C3d, Human, mAb 3 | HM2198 |
| C3g, Human, mAb 9 | HM2199 |
| C4, Mouse, mAb 16D2 | HM1046 |
| C4d, Human, mAb 11D12 | HM2229 |
| C4d, Human, mAb 7H4 | HM2230 |
| C5/C5a, Human, mAb 557 | HM2077 |
| C5/C5a, N-terminus, Human, mAb 561 | HM2076 |
| C5/C5b, Human, mAb 568 | HM2080 |
| C5b-9, Rat, mAb 2A1 | HM3033 |
| C6, Rat, mAb 3G11 | HM3034 |
| Calprotectin, Human, mAb 27E10 | HM2156 |
| Caveolin-1, Rat, mAb 7C8 | HM3014 |
| CD13, Mouse, ER-BMDM1 | HM1083 |
| CD21 (CR2), Human, mAb 21B9 | HM2139 |
| CD34, Mouse, mAb MEC14.7 | HM1015 |
| CD35 (CR1), Human, mAb 31R | HM2107 |
| CD44, Human, mAb NKI-P2 | HM2217 |
| CD73, Human, mAb 4G4 | HM2215 |
| Crry, Rat, mAb TLD-1C11 | HM3032 |
| DPP IV (CD26), Rat, mAb 5E8 | HM3021 |
| Elafin/SKALP, Human, mAb TRAB2F | HM2063 |
| Elafin/SKALP, Human, mAb TRAB2O | HM2062 |
| Elastase, Human, mAb 265-3K1 | HM2174 |
| Endostatin, Human, mAb 1837-46 | HM2188 |
| G-Protein beta 1&2, Human, mAb Raft.1 | HM2166 |
| H-FABP, Human, mAb 66E2 | HM2016 |
| H-FABP, Human, mAb 67D3 | HM2018 |
| JAM-C, Mouse, mAb CRAM-19 H36 | HM1056 |
| L-FABP, Human, mAb K5A6 | HM2051 |
| L-FABP, Human, mAb L2B10 | HM2049 |
| L-Ficolin, Human, mAb GN4 | HM2090 |
| L-Ficolin, Human, mAb GN5 | HM2091 |
| LL37/CAP18, Human, mAb 1-1C12 | HM2071 |
| LL37/CAP18, Human, mAb 3D11 | HM2070 |
| MBL, Human, mAb 3E7 | HM2061 |
| MBL, Human, mAb D8.18 | HM2081 |
| MRP-8, Human, mAb 7C12/4 | HM2175 |
| MRP-14, Human, mAb 1H9 | HM2176 |
| Nectin-2, Mouse, mAb 502-57 | HM1052 |
| Nectin-3, Mouse, mAb 103-A1 | HM1053 |
| NGAL, Human, mAb 697 | HM2193 |
| Nitrotyrosine, mAb HM11 | HM3001 |
| PLVAP, Human, mAb 174/2 | HM2214 |
| Regucalcin, Human, mAb Regucalcin M | HM3018 |
| SAA-1, Human, mAb Reu 86.1 | HM2100 |
| SAA-1, Human, mAb Reu 86.5 | HM2101 |

POLYCLONAL ANTIBODIES

| SPECIFICITY | CAT. # |
|---------------------|--------|
| B-FABP, Human | HP9029 |
| C1q, Rat | HP8021 |
| C3, Mouse | HP8012 |
| C3, Rat | HP8022 |
| C4, Rat | HP8023 |
| C5, Mouse | HP8013 |
| C5L2, Human | HP9036 |
| C5L2, Mouse | HP8015 |
| C5L2, Rat | HP8018 |
| Elafin/SKALP, Human | HP9025 |
| Elastase, Human | HP9027 |
| L-FABP, Human | HP9021 |
| L-FABP, Rat | HP8010 |
| Occludin, Human | HP9047 |



The role of complement in ANCA related necrotizing small vessel vasculitis

C5 depletion with mAb BB5.1

P. Heeringa, PhD, Department of Pathology and Laboratory Medicine, Medical Biology Section, University Medical Center Groningen, Groningen, The Netherlands

Anti-neutrophil cytoplasmic autoantibodies (ANCA) are specific for enzymes in primary granules of neutrophils and peroxidase-positive lysosomes of monocytes. Approximately 90% of ANCA detected in patients with glomerulonephritis and/or vasculitis are specific for proteinase 3 (Pr3) and myeloperoxidase (MPO). ANCA are strongly associated with pauci-immune necrotizing crescentic glomerulonephritis and a spectrum of necrotizing small vessel vasculitis (1).

The pathogenic potential of anti-MPO antibodies has been confirmed in mouse studies (2). These experiments demonstrated that passive administration of murine anti-MPO IgG induces glomerulonephritis similar to that observed in patients.

The effect of pretreatment or intervention with the C5-inhibiting monoclonal antibody BB5.1 was investigated (3). Mice received BB5.1 either 8 hours before or 1 day after disease induction with anti-MPO IgG and lipopolysaccharide. Anti-MPO treated mice developed hematuria, leukocyturia, albuminuria, and crescentic glomerulonephritis. In contrast, BB5.1 pretreatment completely prevented disease development, as evidenced by the absence of urinary abnormalities, a marked reduction in glomerular neutrophil influx and normal renal morphology. Importantly, BB5.1 administration at 1 day after disease induction also resulted in a marked attenuation of urinary abnormalities and a more than 80% reduction in glomerular crescent formation. From these experiments we conclude that inhibition of C5 activation attenuates disease development in the mouse model of anti-MPO IgG-induced glomerulonephritis. Furthermore, these results favor further investigations into the role of complement activation in human MPO-anti-neutrophil cytoplasmic autoantibody-mediated glomerulonephritis, and indicate that inhibition of C5 activation is a potential therapeutic approach in this disease.

- Jennette JC, Xiao H, and Falk RJ: Pathogenesis of vascular inflammation by anti-neutrophil cytoplasmic antibodies. *J Am Soc Nephrol* 2006, 17: 1235-12422
- Xiao H, Heeringa P, Hu P, Liu Z, Zhao M, Aratani Y, Maeda N, Falk RJ, and Jennette JC: Antineutrophil cytoplasmic autoantibodies specific for myeloperoxidase cause glomerulonephritis and vasculitis in mice. *J Clin Invest* 2002, 110: 955-963
- Huugen D, van Esch A, Xiao H, Peutz-Kootstra CJ, Buurman WA, Cohen Tervaert JW, Jennette JC and Heeringa P: Inhibition of complement factor C5 protects against anti-myeloperoxidase antibody-mediated glomerulonephritis in mice. *Kidney Int* 2007, 71: 646-654



New assay for L-Ficolin quantification

Human L-Ficolin ELISA:

L-Ficolin recognizes distinct danger-associated molecular patterns (DAMP) like GlcNAc-structures in LTA and fungal 1,3-β-D-glucan.

L-Ficolin also recognizes N-acetylated carbohydrates and other non-carbohydrate acetylated compounds such as acetylcholine.

Furthermore, L-Ficolin recognizes apoptotic cells and participates in the removal of host cells. L-ficolin circulates in complex with MASP-2 and can activate the lectin pathway.

Low serum levels of L-Ficolin are reported to be associated with recurrent respiratory infections in children.

Interestingly, L-Ficolin has been implicated in the unique immune challenge during pregnancy. In maternal plasma of normal pregnancies a 4- to 5-fold increase in L-Ficolin was detected compared to healthy non-pregnant persons.

However, significantly lower L-Ficolin levels were associated with preeclamptic pregnancies.

Therefore, assessment of L-Ficolin is warranted to study its regulatory role in the innate immune system.

Now available:

| ELISA | |
|---|---------------------|
| SPECIFICITY | CAT. # |
| L-Ficolin, Human | Unique HK336 |
| Functional MBL/MASP-2, Human, (C4 deposition) | Unique HK327 |
| MBL, Human | Unique HK326 |
| MBL, Human, (lectin activity) | HK323 |
| SP-D, Human | HK335 |
| TCC/sC5b-9/MAC, Human | HK328 |

Also available:

- ▶ Antibodies to human MBL, MASP-1, MASP-2, MASP-3, L-Ficolin, H-Ficolin and M-Ficolin.

Soluble MAdCAM-1, an inflammatory biomarker

The unique Human sMAdCAM-1 ELISA (Cat.# HK337) offers a sensitive, non-invasive tool for monitoring chronic inflammatory disease activity in diverse fields of inflammatory research.

The mucosal addressin cell adhesion molecule-1 (MAdCAM-1) is a cell-surface Ig superfamily member of ~60 kDa. MAdCAM-1 is preferentially expressed on the surface of high endothelial venules (HEV) in the gut and associated lymphoid tissue (Peyer's patches). MAdCAM-1 promotes the adhesion of predominantly T and B cells, monocytes/macrophages to the vascular endothelium and is critical for lymphocyte homing to the gut. MAdCAM-1 is upregulated on gut lamina propria in inflammatory bowel disease (IBD), especially Crohn's disease.

Higher expression of MAdCAM-1 is reflected in elevated levels of the circulating soluble form of MAdCAM-1 (sMAdCAM-1). Since MAdCAM-1 is elevated in inflammatory, infectious and malignant diseases, sMAdCAM-1 serves as a perfect non-invasive biomarker for disease activity.


Special features:

- Useful for quantitative measurement of human sMAdCAM-1 in plasma, urine, mother milk and other body fluids
- Standard curve: 0.4 to 100 ng/ml
- Sample volume: less than 25 ul plasma per determination

Hot new products

Anti-TNF α clone 5N, a new strong TNF inhibitory monoclonal antibody

Thirty times stronger TNF α inhibition by monoclonal antibody 5N than by monoclonal T1. Inhibitor properties assessed in Actinomycin D L929 test.



MONOCLONAL ANTIBODIES

| SPECIFICITY | CAT. # |
|-------------------------------|-----------|
| TNF α , Human, mAb, 5N | HM2218 |
| TNF α , Human, mAb, 5N | HM2218-05 |
| TNF α , Human, Protein | HC2040 |