TROUBLESHOOTING IMMUNOFLUORESCENCE



Troubleshooting

This troubleshooting document gives a guideline to the problem, possible cause and suggested solution for problems during the immunofluorescence application:

Problem: \	Weak or no	staining
------------	------------	----------

Problem: Weak or no staining		
The primary antibodies are	Use a new batch of antibodies.	
inactive		
The primary/secondary	Store antibodies according to manufacturer's instructions.	
antibodies may have lost its	Aliquot antibodies into smaller volumes and store in freezer (-20 °C to -	
activity due to improper	70 °C), avoid repeated freeze and thaw cycles.	
storage		
Improper dilution has taken	Run positive controls to make sure that the primary/secondary	
place, or extensive	antibodies are working properly.	
freezing/thawing		
The protein of interest is not	Run an appropriate positive control.	
present in the sample		
The protein of interest is not	Maximize the signal with an amplification step.	
abundantly present in the		
sample		
The antibody concentration is	Increase the concentration of primary/secondary antibodies. Or	
too low	determine the optimal dilution, by running a serial dilution test, to find the	
	best signal to noise ratio.	
The antibody incubation time	Increase antibody incubation time, incubate longer (e.g. overnight) at	
is inadequate	4℃.	
The fixation is inadequate	Try different fixatives or increase the duration of post-fixation.	
The epitope that the antibody	Reduce the duration of post-fixation.	
recognizes may be modified	If the sample has already been overfixed, perform an appropriate or	
by fixation procedures, which	recommended antigen retrieval procedure.	
use formalin and		
paraformaldehyde fixatives		
The antibody may not be	Test the antibody in a native (non-denatured) form on a Western Blotting	
suitable for IF procedures	to ensure it is not damaged.	
which reveal the protein in its		
native (3D) form		
The primary antibody and the	Use secondary antibody that will interact with primary antibody,	
secondary antibody are not	antibodies that are raised against the species in which the primary was	
compatible	raised. For example, if primary antibodies are raised from rabbits, use	
	anti-rabbit secondary antibodies.	
The secondary antibody is	Use a new batch of antibodies.	
inactive	Harana was bakaba at was was ta	
The enzyme substrate system	Use a new batch of reagents.	
is defect or incompatible	La conserva de la contrada del contrada de la contrada de la contrada del contrada de la contrada del contrada de la contrada de la contrada de la contrada del contrada de la contrada del contrada del contrada de la contrada de la contrada del contrada del contrada de la contrada del contrada	
The substrate incubation time	Increase the substrate incubation time.	
is inadequate	Add a payman hilising against to the positional will store by the contribution	
The protein is located in the	Add a permeabilizing agent to the antibody dilution buffer and blocking buffer.	
nucleus (nuclear protein) and	buller.	
the antibody cannot penetrate		
the nucleus	Change correct mounting medium	
Incorrect mounting medium	Choose correct mounting medium.	
The reagents are applied in	Check the protocol used.	
wrong order or steps omitted	Alverse many and the fit research and the distance of the second and the second a	
The fluorescent antibody was	Always prevent the fluorescent-antibody from exposure to light.	
not stored in the dark.		

Problem: Non-specific staining

The concentration of primary/ secondary antibodies is too high	Reduce the concentration of primary/secondary antibodies. Or determine the optimal dilution, by running a serial dilution test, to find the best signal to noise ratio.
The incubation time is too	<u> </u>
long	
The incubation temperature is too high	Reduce the incubation temperature.
The substrate incubation time	Reduce the substrate incubation time.

All Hycult Biotech products are subject to strict quality control procedures.

Copyright® Hycult Biotech. All rights reserved.

The information on this sheet should neither be considered comprehensive or definitive.

TROUBLESHOOTING IMMUNOFLUORESCENCE



is too long	
The primary antibody is raised	Use secondary antibody that will interact with primary antibody,
against the same species as	antibodies that are raised against the species in which the primary was
the sample stained (e.g.	raised. For example, if primary antibodies are raised from rabbits, use
mouse primary antibody	anti-rabbit secondary antibodies.
tested on mouse sample).	Treat sample with e.g. Mouse-On-Mouse blocking reagent prior to the
When the secondary antibody	primary antibody incubation.
is applied it binds to all the	
sample as it is raised against	
that species	
The samples have dried out	Avoid drying out of the samples.

Problem: High Background

Problem: High Background	
The sample is not washed	Wash the sample at least 3 times in PBS between all steps.
enough, the fixative is still	
present.	
The sample contains	Block endogenous enzyme activities using 3% H ₂ O ₂ in methanol or
endogenous enzyme such as	levamisole solution before the incubation of the primary antibodies.
peroxidase or alkaline	
phosphatise	
The sample contains	Block endogenous biotin activity using avidin/biotin blocking reagent
endogenous biotin activity	prior to incubation of primary antibodies.
Blocking of non specific	Increase the blocking incubation period and consider changing blocking
binding might be absent or	agent. 10% normal serum 1hr for sections or 1-5% BSA for 30 min for
insufficient	cells in culture.
The primary antibodies bind	Non-specific binding may be reduced by using higher dilution of the
non-specifically to the sample	primary antibodies.
Antibody concentration was	Titrate the antibody to the optimal concentration, incubate for longer but
too high	in more dilute antibody (a slow but targeted binding is best).
The incubation temperature is	Incubate sections or cells at 4° C.
too high	
The secondary antibodies	Run a secondary control without primary antibody.
bind non-specifically to the	Treat sample with normal serum from the same species as secondary
sample	antibodies, or use pre-adsorbed secondary antibody.
The chromogen reacts with	Use Tris buffer to wash samples.
the PBS present in the cells	
Inadequate fixation causes	Increase the duration of post-fixation.
diffusion of the antigen	
Pemeabilization has damaged	Use buffers without permeabilizing agent.
the membrane and removed	
the membrane protein.	
The primary antibody is raised	Use secondary antibody that will interact with primary antibody,
against the same species as	antibodies that are raised against the species in which the primary was
the sample stained (e.g Mouse	raised. For example, if primary antibodies are raised from rabbits, use
primary antibody tested on	anti-rabbit secondary antibodies.
mouse sample). When the	Treat sample with e.g. Mouse-On-Mouse blocking reagent prior to the
secondary antibody is applied	primary antibody incubation
it binds to all the sample as it	
is raised against that species.	
The sections have dried out	Avoid drying out of the sections.

Helpful links / references

www.protocol-online.org www.ihcworld.com