

Quantitative measurement of methylglyoxal (MGO), a potent and damaging glycating metabolite, by ELISA

Jan van Groningen, Geurt Schilders, Helma W.H. Rutjes, Jan van Binsbergen

Introduction

Methylglyoxal (MGO) is a highly potent glycating metabolite which is produced during glycose, amino acid as well as fatty acid metabolism. MGO forms covalent adducts with the side chains of lysine-, cysteine- and arginine-residues, which is classified as an early glycation process and occurs in all tissues and body fluids. Subsequent, later-stage reactions result in the formation of advanced glycation end products termed AGEs.

The damaging effect of MGO consist of the covalent cross-linking of proteins, or when the modified amino acid is located in the active site of the enzyme or results in a structural rearrangement of the modified protein. In order to counteract the damaging effect of MGO our body has evolved several non-enzymatic as well as non-enzymatic defensis such as the glyoxalase system. In C. elegans silencing of the glyoxalase system resulted in a decreased lifespan by 40%. In human, increased production and accumulation of methylglyoxal (MGO) are hallmarks of aging and a pathological conditions number of such as neurodegenerative disorders and diabetes.

In order to quantitatively measure MGO an ELISA has been developed using MGO-specific antibodies



Features of the assay

- Assay time: 3.5 hours
- · Species independent
- · Sample types tested: plasma, serum, urine and faeces
- HIT503 standard consists of MGO labelled HSA (µg/ml)
- Detection limit: 7.8µg/ml HSA-MGO equivalents.
- Dynamic range: 7.8 500µg/ml HSA-MGO
- equivalents.
- Sample volume: 125µl/well
- Non-expensive

Skeletal formula of MGO



Measurement of MGO in human EDTA plasma

A panel of 20 plasma samples from randomly selected healthy donors (males and females) was tested for MGO levels
The average level of MGO was 456µg/ml HSA-MGO equivalents, with a range of 358-636µg/ml HSA-MGO eq.



Detection of MGO with HIT503 assay is species independent



MGO levels were measured in sera of various species using the HIT503 assay The calculated MGO concentrations are depicted in Table 4.

	HSA-MGO eq (µg/ml)	SD	CV (%)
dog	504	42	8
horse	589	24	4
goat	218	16	8
swine	537	19	3
ray	147	14	10
mouse	145	4	2
rabbit	176	11	6
bovine	533	80	15

www.hycultbiotech.com info@hycultbiotech.com

HIT503 assay characteristics Typical standard curve of inhibitory HIT503 assay Freeze-thaw stability

plasma

3.0		
2.5		
o 2.0		
¥ 1.5	\sim	
° 1.0	· · · · · · · · · · · · · · · · · · ·	\sim
0.5		\sim
0.0		
10	100	100
	HSA-MGO eq. (µg/ml)	

Recovery

- 100 and 500µg of HSA-MGO eq. were spiked into citrate plasma samples (see Table 1)
- Percentage recovery was calculated as the percentage of measured/expected MGO
- Recovery in both plasma samples is between 80-120% Table 1. Recovery of CML in human citrate plasma

matrix	spike (µg/ml)	% recovery	
Citrate plasma 1	500	107	
	100	105	
Citrate plasma 2	500	105	
	100	109	

Table 3. Bench top stability of standard and plasma samples 10 min RT 2 hrs RT 24 hrs 4°C 100% 101% 95% plasma 98% standard 100% 94% 92% 83%

Stability in plasma sample and standard is between 80-

Stability of MGO in human citrate plasma and standard after several freeze-thaw cycles (Table 2).

Stability in plasma is between 80-120%. Stability of standard drops after repeated freeze-thaw cycles. Table 2. Freeze-thaw stability of standard and plasma sam

100% 100%

standard 100% 98% *Number of freeze-thaw cycles

Bench top stability

120%

1x 2x 3x 4x 5x*

93% 94% 94%

84%

Recovery of MGO in human citrate plasma sample and standard after incubation for 10 min. at RT or 2 or 24hrs at 4°C as compared to freshly prepared samples.

81%

62%

Detection of MGO with HIT503 assay in different matrices

MGO levels were measured in plasma, urine and faeces samples using the HIT503 assay Linearity of dilution was determined by serially diluting samples. The "measured" conc "expected" concentrations. A correlation of at least 0.97 was found in the matrices tested. ns are plotted against the



Hycult Biotech (HQ) PO Box 30 5400 AA Uden The Netherlands Tel: +31-413-251335

Hycult Biotech Inc. 600 West Germantown Pike Suite 400 Plymouth Meeting, PA 19462. USA Tel: 1-855-2-HYCULT (855-249-2858) 610-260-1491 (local)