

BIOTINYLATED RAT IgG_{2b} ISOTYPE CONTROL

Catalog nr	HI3008 (lot number and expiry date are indicated on the label)	
Description	A rat IgG _{2b} negative control for the measurement of non-specific binding of rat monoclonal antibodies of isotype IgG _{2b} . The control antibody is suitable for use as a negative control with human, mouse and canine tissues.	
Species	Rat IgG _{2b}	
Formulation	0.5 ml (100 µg/ml) 0.2 µm filtered biotinylated antibody solution in PBS, containing 0.1% bovine serum albumine and 0.02% sodium azide.	
Application	The biotinylated antibody can be used for flow cytometry.	
Use	For flow cytometry dilutions to be used depend on detection system applied. It is recommended that the user titrates the antibody for use in their system to a concentration equivalent to that of their test reagents.	
Storage and stability	Product should be stored undiluted at 4°C. Under recommended storage conditions, product is stable for one year.	
Precautions	For research use only. Not for use in or on humans or animals or for diagnostics. It is the responsibility of the user to comply with all local/state and federal rules in the use of this product. Hbt is not responsible for any patent infringements that might result from the use or derivation of this product.	
Also available	HI3007	Rat IgG _{2b} isotype control
	HI3009	FITC conjugated Rat IgG _{2b} isotype control
	HI3002	Biotinylated Rat IgG ₁ isotype control
	HI3005	Biotinylated Rat IgG _{2a} isotype control