

CERTIFICATE OF ANALYSIS – TECHNICAL DATA SHEET

Product name Mannose receptor, Mouse, clone MR5D3

Catalog number HM1049-100UG

Lot number **Expiry date**

100 μg Volume 1 ml **Amount**

Formulation 0.2 µm filtered in PBS+0.1%BSA Concentration 100 μg/ml

Host Species Rat IgG2a Conjugate None

Endotoxin <24 EU/mg Purification Protein G

4°C Storage

Application notes

	IHC-F	IHC-P	IF	FC	FS	IA	IP	W
Reference #	2,3,4		6,9	2,6,8-10	6	1,5,7	2	1,2
Yes	•		•	•	•	•	•	•
No								

N.D.= Not Determined; IHC = Immuno histochemistry; F = Frozen sections; P = Paraffin sections; IF = Immuno Fluorescence; FC = Flow Cytometry; FS = Functional Studies; IA = Immuno Assays; IP = Immuno Precipitation; W = Western blot

Dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:50.

- FC: 2% paraformaldehyde fixed cells were labeled with 10µg/ml of mAb.
- IHC-F:Sections fixed in ethanol or 2% paraformaldehyde were permeabilized with PBS/0.1% triton X100. After quenching and blocking mAb was added (10μg/ml) and incubated for 1h
- IA: mAb was used as detector at 10µg/ml for 2h.
- IP:MR was immunoprecipiated using antibody MR5D3 (10µg/ml) and GammaBind Plus sepharose (Ref. 2)
- W: Cleared lysates were run on 6% SDS-PAGE under nonreducing conditions. After transfer blot was stained with mAb (2μg/ml). FS: Antibody MR5D3 was used as a MR ligand to target MR+ dendtritic cells resulting in anti-rat IgG production (Ref. 6).

General Information

Description

The monoclonal antibody MR5D3 recognizes the Mannose Receptor (MR, also known as CD206), a member of the vertebrate C-type lectin family. MR is a pattern recognition receptor that is involved in both innate and adaptive immunity. The 175 kDa transmembrane protein consists of 5 domains: an amino-terminal cysteine-rich region, a fibronectin type II repeat, a series of eight tandem lectin-like carbohydrate recognition domains (responsible for the recognition of mannose and fucose), a transmembrane domain and an intracellular carboxy-terminal tail. The structure is shared by the family of multi lectin mannose receptors: the phospholipase A2-receptor, DEC 205 and the novel Ctype lectin receptor (mannose receptor X). The MR recognises a wide range of gram positive and gram negative bacteria, yeasts, parasites and mycobacteria. Moreover, the MR has also been shown to bind and internalize tissuetype plasminogen activator.

MR's are present on monocytes and dendritic cells (DC) and are presumed to play a role in innate and adaptive immunity, the latter via processing by DC. The expression of MR as observed in immunohistology is present on tissue macrophages, dendritic cells, a subpopulation of endothelial cells, Kupffer cells and sperm cells. The expression of MR on monocytes increases during culture and can be enhanced by cytokines such as IL-4. Labelling of MR expressing monocytes/macrophages increases with prolonged incubation time probably due to internalization of the MR-antibodycomplex. Detection of the MR with anti-MR monoclonal antibody MR5D3 can substitute staining for mannose containing probes as labeled mannosylated BSA, a technique which is more cumbersome and less specific.

Chimaeric CRD4-7-Fc protein

Aliases CD206, MMR, cLEC13D, Macrophage mannose receptor 1

References

Immunogen

- Zamze, S et al; Recognition of bacterial capsular polysaccharides and lipopolysaccharides by the macrophage mannose receptor. J Biol Chem 2002, 277: 41613
- 2. Martinez-Pomares, L et al. Analysis of mannose receptor regulation by IL4, IL10 an proteolytic processing using novel monoclonal antibodies. J. of Leuk Biol 2003, 73:604
- 3. Hassan, M et al. The schistosoa mansoni hepatic egg granuloma provides a favourable microenvironment for sustained growth of leishmania donovani. Am J Path 2006, 169:943
- Duspasquir, M et al. The dermal microenvironment induces the expression of the alternative activation marker CD301/mGL in mononuclear phagocytes, independent of IL4/IL13 signalling. J Leuk Biol 2006, 80:838

Version: 11-2019

- Boskovic, J et al. Structural model for the mannose receptor family uncovered by electron microscopy of endo 180 and the mannose receptor. J Biol Chem 2006, 281:8780
- McKenzie, E et al. Mannose receptor expression and function define a new population of murine dendritic cells. J immunol 2007, 178:4975
 Miller, J et al. The mannose receptor mediates dengue virus infection of microphages. Plos pathogens 2008,
- 7. 4:0001
- Segura, E et al. Different cross-presentation pathways in steady-state and inflammatory dendritic cells. PNAS
- Hardison, S et al. Interleukin-17 is not required for classical macrophage activation in a pulmonary mouse model of Cryptococcus neoformans infection. Infect&immune 2010, 78:5341
- Geier, H et al. Phagocytic receptors dictate phagosomal escape and intracellular proliferation of francisella tularensis. Infect Immun 2011, march 21

Storage&stability

Product should be stored at 4°C. Under recommended storage conditions, product is stable for at least one year.

Precautions

For research use only. Not for use in or on humans or animals or for diagnostics. It is the responsibility of the user to comply with all local/state and federal rules in the use of this product. Hycult Biotech is not responsible for any patent infringements that might result from the use or derivation of this product.

We hereby certify that the above-stated information is correct and that this product has been successfully tested by the Quality Control Department. This product was released for sale according to the existing specifications. This document has been produced electronically and is valid without a signature.

Approved by Manager of QC Brenda Teunissen

Date 12/11/2019

Do you have any questions or comments regarding this product? Please contact us via support@hycultbiotech.com.

Version: 11-2019