

CERTIFICATE OF ANALYSIS – TECHNICAL DATA SHEET

Product name	IFN-Gamma, Rat, Polyclonal		
Catalog number	HP8007		
Lot number	-	Expiry date	-
Volume	1 ml	Amount	100 µg
Formulation	0.2 µm filtered in PBS+0.1%BSA	Concentration	100 µg/ml
Host species	Rabbit IgG	Conjugate	None
Endotoxin	<24 EU/mg	Purification	Protein A
Storage	4°C		

Application notes

	IHC-F	IHC-P	IF	FC	FS	IA	IP	W
Reference #								
Yes	•				•	•		•
No								
N.D.		•	•	•			•	

N.D.= Not Determined; IHC = Immuno histochemistry; F = Frozen sections; P = Paraffin sections; IF = Immuno Fluorescence; FC = Flow Cytometry; FS = Functional Studies; IA = Immuno Assays; IP = Immuno Precipitation; W = Western blot

Dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:10.

- FS: For neutralization of biological activity dilutions have to be made according to the amounts IFN-gamma to be inactivated.

General Information

Description	The antibody reacts with rat Interferon gamma (IFN-gamma) of both natural and recombinant origin. IFN-gamma is a pluripotent cytokine with important pro-inflammatory functions. The antibody inhibits the biological activity of natural and recombinant rat IFN-gamma. The antigen specificity was further assessed by ELISA and Western blotting. No cross-reactivities with other cytokines have been detected.
Aliases	Interferon gamma
Storage&stability	Product should be stored at 4°C. Under recommended storage conditions, product is stable for at least one year.
Precautions	For research use only. Not for use in or on humans or animals or for diagnostics. It is the responsibility of the user to comply with all local/state and federal rules in the use of this product. Hycult Biotech is not responsible for any patent infringements that might result from the use or derivation of this product.

We hereby certify that the above-stated information is correct and that this product has been successfully tested by the Quality Control Department. This product was released for sale according to the existing specifications. This document has been produced electronically and is valid without a signature.

Approved by Manager of QC
Robbert Zwinkels

Date
08/05/2018

Do you have any questions or comments regarding this product? Please contact us via support@hycultbiotech.com.