

**CERTIFICATE OF ANALYSIS – TECHNICAL DATA SHEET**

<b>Product name</b>	LBP, Human, pAb	<b>Expiry date</b>	-
<b>Catalog number</b>	HP9023		
<b>Lot number</b>	-	<b>Amount</b>	100 µg
<b>Volume</b>	1 ml	<b>Concentration</b>	100 µg/ml
<b>Formulation</b>	0.2 µm filtered in PBS+0.1%BSA+0.02%NaN3	<b>Conjugate</b>	None
<b>Host Species</b>	Rabbit IgG	<b>Purification</b>	Protein A
<b>Endotoxin</b>	N.A.		
<b>Storage</b>	4°C		

**Application notes**

	IHC-F	IHC-P	IF	FC	FS	IA	IP	W
Reference #								
Yes						•	•	•
No								
N.D.	•	•	•	•	•			

N.D.= Not Determined; IHC = Immuno histochemistry; F = Frozen sections; P = Paraffin sections; IF = Immuno Fluorescence; FC = Flow Cytometry; FS = Functional Studies; IA = Immuno Assays; IP = Immuno Precipitation; W = Western blot

Dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:10.

- IA: the antibody is useful to detect LBP in ELISA both as coating and as detector.

**General Information**

<b>Description</b>	The antibody reacts with human natural and recombinant lipopolysaccharide binding protein (LBP). LBP is an approximately 60 kDa acute phase protein that is produced by hepatocytes. This protein strongly binds to LPS and has been shown to play an important role in the handling of LPS by the host. A number of functions of LBP have been reported. First, LBP transfers LPS to the LPS receptor CD14 on mononuclear phagocytes, leading to an 100-1,000-fold increased sensitivity of the cells to LPS. Furthermore, LBP can enhance the response of CD14 negative cells by acceleration of LPS binding to soluble CD14, a complex that stimulates these cells. Next, LBP transfers LPS into High Density Lipoprotein (HDL), which effectively neutralizes its biological potency. LBP was demonstrated to protect mice from septic shock caused by LPS or gram negative bacteria.
<b>Aliases</b>	Lipopolysaccharide Binding Protein
<b>Storage&amp;stability</b>	Product should be stored at 4°C. Under recommended storage conditions, product is stable for at least one year.
<b>Precautions</b>	For research use only. Not for use in or on humans or animals or for diagnostics. It is the responsibility of the user to comply with all local/state and federal rules in the use of this product. Hycult Biotech is not responsible for any patent infringements that might result from the use or derivation of this product.

We hereby certify that the above-stated information is correct and that this product has been successfully tested by the Quality Control Department. This product was released for sale according to the existing specifications. This document has been produced electronically and is valid without a signature.

Approved by Manager of QC  
 Robbert Zwinkels

Date  
 13/03/2018

Do you have any questions or comments regarding this product? Please contact us via [support@hycultbiotech.com](mailto:support@hycultbiotech.com).